# VERMEIREN

## Kim

USER MANUAL  
MANUEL D'UTILISATION  
GEBRUIKSAANWIJZING  
GEBRAUCHSANWEISUNG  
MANUALE DI ISTRUZIONI  
MANUAL DE INSTRUCCIONES  
INSTRUKCJA OBSŁUGI  
NÁVOD K OBSLUZE

EN

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*Instruction manual*

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*Manual de instrucciones*

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This page is intentionally left blank

**EN**

**Contents**

1.	Technical details	2
2.	Components	2
3.	Explanation of symbols	2
4.	Check after delivery	2
5.	Intended use	2
6.	Assembly / Dismantling	2
7.	Folding / Unfolding pedal exerciser	3
8.	Mounting calorie meter	3
9.	Working calorie meter	3
10.	Exercises of the legs	4
11.	Resistance adjustment	4
12.	Safety instructions	4
13.	Maintenance, inspection, cleaning	4

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### **Pedal Exerciser Kim**

#### 1. Technical details

**⚠ WARNING: Risk of unsafe settings - Use only the settings described in this manual.**

Brand	Vermeiren	<b>CE</b>
Type	Pedal exerciser	
Model	Kim	
Width	400 mm	
Length	500 mm	
Height	300 mm	
Folded width	400 mm	
Folded length	300 mm	
Folded height (calorie meter, adjustable star button removed, legs folded)	170 mm	
Weight	2,35 kg	
Chargeable weight	max. <b>50 kg</b>	
Storage temperature	+5°C - +41°C	
We reserve the right to introduce technical changes. Measurement tolerance +/- 15 mm / 1,5 kg		

First of all we wish to thank you for the trust you placed in us by selecting a **VERMEIREN** product.

Before using the pedal exerciser, read the instruction manual carefully: it will familiarise you with this product.

If you have any further questions, please consult your specialist dealer. To find a service facility or specialist dealer near you, contact the nearest Vermeiren facility.

Visually impaired people can contact the dealer for the instructions for use.

#### 2. Components

1. Main frame
2. Pedals
3. Adjustable star button
4. Locking pin folding
5. Calorie meter
6. Identification plate

#### 3. Explanation of symbols

Maximum mass

CE conformity

Type designation

#### 4. Check after delivery

Unpack your product and verify if the delivery is complete. Following items should be included:

- Pedal exerciser frame with 2 pedals
- Adjustable star button (resistance)
- Calorie meter (without battery)
- Manual

Verify your product for transport damage. If you find any damages after delivery, proceed as follows:

- contact the transporter
- have a list made of any problems
- contact your supplier

#### 5. Intended use

The pedal exerciser shall be used by 1 person.  
The pedal exerciser shall be used indoors.

The construction of the pedal exerciser allows full use by elderly people or during revalidation in home environment. The pedal exerciser can increase the flexibility by exercising the range of motion of the joints of the lower body (exercise of the legs).

The pedal exerciser shall not be used as a toy for children.

The pedal exerciser should only be used on horizontal, dry and non-smooth surfaces.

The manufacturer is not liable for damage caused by the lack of or improper service, or as a result of not following instructions from this manual.

#### 6. Assembly / Dismantling

##### 6.1. Assembly pedal exerciser

- Unfold the pedal exerciser.
- Mount the adjustable star button on the pedal exerciser.
- Mount the calorie meter on the pedal exerciser.


##### 6.2. Dismantling pedal exerciser

- Remove the adjustable star button from the pedal exerciser.
- Remove the calorie meter from the pedal exerciser.
- Fold the pedal exerciser.

#### 7. Folding / Unfolding pedal exerciser

**⚠ CAUTION: Risk of clamping – Keep fingers, cloths, jewellery or items away from moving parts of the pedal exerciser (folding legs and upper frame when folding the pedal exerciser).**

To transport the pedal exerciser in the car, or for storage of the pedal exerciser, you can fold the pedal exerciser.


### 7.1. Unfolding pedal exerciser

- Fold the legs (2) outwards until the locking pin (1) clicks automatically in the hole.
- Check that the legs (2) are firmly fixated and cannot move.

### 7.2. Folding pedal exerciser

- Press the locking pin (1).
- Fold the legs inwards (2).

## 8. Mounting calorie meter

- Place the star button (3) through the attachment plate of the calorie meter (4).
- Mount the star button and calorie meter on the pedal exerciser.
- Attach the clasp (5) of the calorie meter to the frame tube (6).

## 9. Working of the calorie meter


- Press the red selection button to activate the calorie meter, or start pedalling with your pedal exerciser. When the magnet passes by the calorie meter, the calorie meter activates automatically.
- The calorie meter display turns off automatically after 4 to 5 minutes.
- When pedalling with the pedal exerciser, the STOP in the bottom right corner disappears, and all the functions start calculating.
- If you stop pedalling, the STOP will reappear in the bottom right corner after 4 seconds, all functions stop calculating.
- The calorie meter displays 4 different functions.
  - CNT:** Gives the amount of revolutions of the pedals (counts how many times the magnet on the pedal passes the calorie meter) from this exercise session.
  - TOT.CNT:** Gives the total amount of revolutions of the pedals from different exercise sessions.
  - CAL:** Gives the total amount of calories consumed during the exercise session.
  - TIME:** Gives the total time of the exercise session.

Press the red selection button to go to the next function (TOT.CNT, CAL, TIME). The display always shows the CNT function.

Press the red selection button for 4 seconds to reset the functions (CNT, CAL, TIME). The TOT.CNT is not reset. When you have reset the data, a next exercise session shall be started.

- **SCAN:** Press the red selection button to go to SCAN. When the arrow points to SCAN, you have selected the mode SCAN. Now the display automatically goes through all the functions (TIME, TOT.CNT, CAL), and shows them consecutively.
- When you remove the battery, all the data (TIME, TOT.CNT, CNT, CAL) shall be lost. If you see an improper display, please replace the battery to have a good result. Use a 1,5 V battery as power supply.

## 10. Exercises of the legs

**⚠ CAUTION: Risk of injury – Always sit in a chair when performing the leg exercise. Do not stand on the pedal exerciser.**

**⚠ CAUTION: Risk of injury - Make sure that the chair does not move or rotate. If the chair has wheels, be sure the brakes are engaged.**

**⚠ CAUTION: Risk of injury - Make sure that you wear good shoes that do not slip off the pedals. (Rubber soled shoes)**

**⚠ CAUTION: Risk of burns – Do not take the pedal exerciser by the pedals tube after training, this could be warm caused by the friction of moving parts.**

- Place pedal exerciser on a horizontal, non-smooth surface directly in front of the chair.
- Sit comfortably in a non-movable chair.
- Place your feet in the pedals.
- Now you can pedal with the pedal exerciser.

## 11. Resistance adjustment

Turn the adjustable star button to adjust the resistance for pedalling.

- Turn clockwise (direction +): Resistance increases -- heavier pedalling.
- Turn counter clockwise (direction -): Resistance decreases -- looser pedalling.

## 12. Safety instructions

- Consult your specialist (doctor) before using the pedal exerciser.

- Stop the exercises and consult your specialist immediately if you feel dizziness, severe muscle/joint pain or chest pain.
- Warm-up and stretch about 5 minutes before beginning your exercises.
- Exercises shall be in a controlled manner and beginning slowly.
- Do not exhaust yourself.
- **Do not exceed the chargeable weight. The chargeable weight for your pedal exerciser is 50 kg.**
- Use the pedal exerciser on horizontal, dry, non-smooth and indoor surfaces.
- Keep a minimum safety area clearance of 1 meter around the pedal exerciser. Make sure that this area remains clear during use of the pedal exerciser.
- Keep away from children. The pedal exerciser shall not be used as a toy.
- When it is used by or near children, invalids or disabled persons, close supervision is necessary.
- Risk of burns – Be careful when using in hot or cold environments (sunshine, extreme cold, saunas, etc.) and when touching.
- Do not place fingers between components of the pedal exerciser when folding.
- Keep fingers, clothes, jewellery or other items away from moving parts of the pedal exerciser (rotating pedals).
- Wear comfortable and suitable clothing when using the pedal exerciser.
- Be sure to wear good shoes that do not slip off the pedals.
- Make sure that the chair in which you are sitting does not move or rotate.
- Do not stand on the pedal exerciser.
- Do not use the pedal exerciser when it is not working properly.
- Take note of the instructions for care and service. The manufacturer is not liable for damage caused by improper servicing/care.

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## Service registration form

This product (name): .....

was inspected (I), serviced (S), repaired (R) or disinfected (D):

By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:
By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:
By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:
By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:
By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:	By (stamp):  Kind of work: I / S / R / D Date:

Vermeiren GROUP  
Vermeirenplein 1 / 15  
B-2920 Kalmthout  
BE

website: [www.vermeiren.com](http://www.vermeiren.com)

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### All rights reserved, including translation.