



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

<h3 style="text-align: center;">Differential for 2217 with Free Wheel</h3>				---
				Plan No. EP00223/G
Scale: 1:5	Drafted	Verified	Confirmed	
Date	01.12.2016	28.11.2025	28.11.2025	
Initials	KR	PM	AZ	

Material			
No.	Plan	Article	Quantity
1	YP04707	2886600	2
2	---	---	-
3	X002282	1992854	6
4	6x6x50 DIN 6885A	GP19000	2
5	XP02369	GP17400	1
6	XP02368	GP17300	1
7	XP01189	2886004	2
8	XP01195	2992221	2
9	XP01197	2992222	2
10	XP01202	2992224	1
11	XP03236	2631200	1
12	XP01283	2992321	1
13	XP01319 (23T)	GP06200	1
	XP03162 (22T)	GP25300	1
14	D10,5 DIN 125	1001558	3
15	M10 DIN 934	1001672	4
16	Chain	2993196	1,4 [m]
17	D15 DIN 6799	1003049	4
18	(23T)	1008502	1
	XP03160 (22T)	2996728	1
19	EF00100 Front Wheel 20"	2991678	1
20	EF00101 Wheel 20"	2994423	2
21	M6 DIN 127B	1004414	4
22*	---	1902733	3
23	M8x25	1990501	2
24	M5x16	1006112	6
25	---	1994568	6
26	D8,4 DIN 9021	1996574	2
27	M6x90 DIN 912	1911103 (01,02,03)	4
	M6x110 DIN 912	1911104 (06)	4
	M6x100 DIN 912	1998869 (07)	4

Material			
No.	Plan	Article	Quantity
28	M10 DIN 493B	1903216	1
29	XP02365	GP17200	1
30	XP01259	2992306 (01)	1
	XP01671	2993508 (02)	1
	XP01672	2993509 (03)	1
	XP01673	2993510 (06)	1
	XP01674	2993511 (07)	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Plan No.
EP00223/G


Differential for 2217 with Free Wheel			
Scale: 1:5	Drafted	Verified	Confirmed
Date	01.12.2016	28.11.2025	28.11.2025
Initials	KR	PM	AZ

* not depicted