




Sheet 1/2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

VERSO SE53 BI5 EXPLOSION

Scale:	2:5	Drafted	Verified	Confirmed
Date	4/03/2025	9/04/2025	9/04/2025	
Initials	SD	PM	MK	

Plan No.

E00100260/D.1



Material			
No.	Plan	Article	Qty
1	P00052550	2085754	1
2	P00096804	2093354	3
3	P00096805	2093454 (R)	1
	P00096980	2093554 (L)	1
4	P00096810	2093654	1
5	P00096975	2093754 (R)	1
	P00096976	2093854 (L)	1
10	M6 DIN125	1000973	8
11	20x20	1001612	5
12	M10 DIN934	1001672	1
13	---	1919578	2
14	M10x19	1991432	1
15	M6 DIN985	2994063	4
16	M6x40 ISO7380	2994066	4
17	M6 DIN125	2994414	8
18	M5x16 DIN912	1915415	2
19	M6x10	2999405	1
20	M5x25 DIN912	2994906	2
21	M5 DIN127	2995288	2
22	---	1990752	1

| Electrical components are presented in 2998566

Sheet 2/2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

VERSO SE53 BI5 EXPLOSION			
Scale: 2:5	Drafted	Verified	Confirmed
Date	4/03/2025	9/04/2025	9/04/2025
Initials	SD	PM	MK

Plan No.
E00100260/D.1

