



SE53 WIDE ARMREST / TURIOS

Plan No.

EP00489/ I



Material				Material				Bill of Material 1913465/1913464 for EP00489			
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity	No.	Plan	Article	Quantity
1	-	2606454	1	28*	M4X12	2994126	2	11a	X005567	1912286 (R)	1
2	Y003626	1476554	1	29	Linx DLX REM 050	1909020	1		X005566	1912285 (L)	1
3	Y009175	1641654	3		Linx DLX REM 060	1910700	1	11b	X005517	2014999	1
4	Y009176	1641754	1		Linx SE19 SE39	1912296	1	11c	X005561	2021699	1
5	YP08065	2605954	2		Linx SE39	1912301	1	11d	M6x16	1008459	1
6A	YP08058	2605454 (L)	1		Linx SE19 ACT5	1912305	1	11e	Nut M6	1008492	1
	YP08059	2605554 (R)	1		Linx LCD +JOYST	1913523	1	11f	Rivet Ø4	1992805	2
6B*	YP08060	2605654 (L)	1	30	M5x16	1907493	2	11g	M5x16	3991136	2
	YP08061	2605754 (R)	1	31*	Light controller	1921014	1	11h	Washer M5	1906561	2
7	-	-	-	32	Spring plunger	1991432	1	11i	Sticker	T001002	1
8	YP09265	2642654	1	33	-	-	-				
9	-	-	-	34a	Armrest wide Soft Gray	1999993 (R)	1				
10	X005565	1912287	1		Armrest wide Soft Gray	1999992 (L)	1				
11	Armrest right	1913465	1		Armrest wide Soft Black	2998254 (R)	1				
	Armrest left	1913464	1		Armrest wide Soft Black	2998253 (L)	1				
12	XP00227	2990147	2/4**	34b	L30	1998790 (R)	1				
13	-	-	-		L30	1998789 (L)	1				
14	Washer M6	1000973	8	35*	Washer M4 DIN127B	1998515	2				
15	M5x12	1001489	2	36	M6x45	2994058	2				
16	Tube Cap 20x20	1001612	2	37	Tube Cap 16x16	2994195	2				
17	Nut M10	1001672	1	38	XP03080	2996656	1				
18	Nut M5	1001679	6								
19	Nut M6	1001680	6								
20	M6x40	1003111	2								
21	Star Knob	1003658	1								
22	M5x20	1003907	4								
23	M6x40	1006111	4								
24	-	-	-								
25	M6x16	1008459	2								
26	Nut M6	1008492	2								
27	-	-	-								

EP00489 / I

SE53 WIDE ARMREST / TURIOS

Scale: 1:5

Drafted

Verified

Confirmed

Date

12/02/2020

13/03/2025

13/03/2025

Initials

DR

PM

MK

Plan No.

EP00489 / I



* - used only with SE19

** - for L30