




Verso			
No.	Plan	Description	Quantity
1	E00054056	Side frame BI5	2
	E00057375	Side frame B02/B03	2
Verso 30°			
No.	Plan	Description	Quantity
1	E00080121	Side frame BI5	2
	E00083256	Side frame BI5 Gas Spring	2
	E00080122	Side frame B02/B03	1
	E00083931	Side frame B03 Gas Spring	1
Verso/Verso 30°			
No.	Plan	Description	Quantity
2	E00054636	Rear frame (44-50)	1
	E00082136	Rear frame (39-42)	1
	E00058416	24 INCH Wheel adaptor	2
3	E00053247	Lock mechanism	2
4	E00054057	Front fork	2
5	E00071954	Footrest B06	2
6	E00077326	BK6	1
	E00055449	Cross	1
7	E00055516	Joystick attachment	1
7	E00075959	SE53	1
	E00057691	Upholstery	1
9	E001038	Armrest BI5	2
10	E00093611	SE23 Armrest extension	2
	2998597	Drive train standard	1
11	2998566	Drive train + SE19 lights	1
	1918834	Battery 24V 12Ah	1
11	1919283	Battery 24V 17Ah	1
12*	1919284	Charger	1
13	E00083496	SE19 LIGHTING	1
14*	E00065443	Wheel cover SE46	1
15*	E00078517	SE08 With wrist support	1

* Not depicted

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

5146249

VERSO/VERSO 30° SPECIFICATION

Plan No.

E00057683/D.1

Scale: 1:8

Drafted

Verified

Confirmed

Date

11/04/2023

15/04/2025

15/04/2025

Initials

SD

AS

MK

