



Material					
No.	Article	XS	S	M	L
1	2991014	0,033	0,055	0,068	0,11
2	1009516	5,4	0	0	0
	1009517	0	5,4	5,4	5,4
3	1991236	12	0	0	0
	1009512	0	12	12	12
4	W700813	5	0	0	0
	1009508	0	5	5	5
5	1992533	5	0	0	0
	2991330	0	5	5	5
6	2991886	0,8	1	1,2	1,4

Sizes [cm]					
	A	B	C	LENGHT	Article
XS	2,5	4,4	15	117	2990680
S	3,8	5,7	18	127	2990681
M	3,8	6,4	23	142	2991098
L	3,8	7,6	28	152	2991099

VERMEIREN NV
 Vermeirenplein 1-15
 2920 Kalmthout
 www.vermeiren.com

Made in Poland
NeoFlex
U75Neo.
size M

contains:
 90% Neopren
 8% Poliester
 2% Elastan



Mounting			
No.	Plan	Article	Quantity
1	XP01976	2993641	4
2	XP01980	2993813	4
3	XP01979	2993812	4
4	M6X30 collar bolt	2994294	4
5	YP07178	2760099	4

EP00345 / B

NEO U75 ABDUCTION BELT 4 POINT				Plan No.
				EP00345/ B
Scale: N/A	Drafted	Verified	Confirmed	
Date	31/07/2018	04/03/2019	04/03/2019	
Initials	FD	AS	TK	