



**VERMEIREN
NV**

Vermeirenplein 1-15
2920 Kalmthout
www.vermeiren.com

Made in Poland

NeoFlex

**U173Neo.
size L**

contains:


- 76% Polyethylene
- 15% Polyester
- 7% Polyethersulfone
- 2% Elastane

SGS INSTITUT
FRESENIUS



Material			
No.	Article	S	L
1	2994147	0,15	0,2
2	2991886	3,5	4
3	2994164	2	2,5
4	2991191	2	2
5	2991330	2	2
6	1995604	0,2	0,2
7	1995605	1	1

Sizes [cm]							
	A	B	C	D	E	F	Article
S	41	20	27	6	25	50	2994151
L	51	20	35	6	32	62.5	2994152



NEO U173 TORSO VEST ZIP

Plan No.

EP00311 / A

Scale: N/A	Drafted	Verified	Confirmed
Date	31/07/2018	31/07/2018	31/07/2018
Initials	FD	AS	TK



EP00311 / A