



A (2:7)

Use Loctite 582 for 22 /
Uzyc Loctite 582 dla 22.

Use Loctite 582 for 23/
Uzyc Loctite 582 dla 23.

B (2:7)

Use lubricant between 34,35/
Uzyc smaru pomiedzy 34,35.

Use Loctite 582 for 22 /
Uzyc Loctite 582 dla 22.

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Differential for 2219 with Nexus 3-speed			
Scale:	1:7	Drafted	Verified
Date	26.05.2025	2.12.2025	2.12.2025
Initials	AZ	PM	AZ

Plan No.
EP00169/J

Material			
No.	Plan	Article	Quantity
1	YP04707	2886600	2
2	XP02365	GP17200	1
3	X002282	1992854	6
4	6x6x50 DIN6885A	GP19000	2
5	XP02371	GP17600	1
6	XP02370	GP17500	1
7	XP01176	2992208	1
8	XP01177	2992207	1
9	XP01189	2886004	2
10	XP01202	2992224	1
11	XP03236	2631200	1
12	XP01283	2992321	2
13	D10,5 DIN 125	1001558	3
14	M10 DIN 934	1001672	4
15	M6 DIN 985	1001680	4
16	D15 DIN 6799	1003049	4
17	M4x16	2990185	6
18	XP03143	2996710	1
19	-----	-----	-
20	Front Wheel 24"	2992350	1
21	D160mm	2992382	1
22	M8x25	1990501	2
23	M5x16	1006112	6
24	---	1994568	6
25	D8,4 DIN 9021	1996574	2
26	M10 DIN 439B	1903216	1
27	M4 DIN 985	1001678	6
28	M6x90 DIN 912	1911103 (01,02,03)	4
	M6x110 DIN 912	1911104 (06)	4
	M6x100 DIN 912	1998869 (07)	4
29	EF00105 Wheel 24" Drive	2994425	2
30*	---	---	-
	Vintage Tyre	1907965	3
31**	---	1902733	3

Material			
No.	Plan	Article	Quantity
32	Chain	2993196	1,66 [m]
33	XP01195	2992221	2
34	XP01197	2992222	2
35	---	G040600	1
36	XP01259	2992306 (01)	1
	XP01671	2993508 (02)	1
	XP01672	2993509 (03)	1
	XP01673	2993510 (06)	1
	XP01674	2993511 (07)	1

Nexus 3 subparts - 2996710 (2992142***)			
No.	Plan	Article	Description
18a	YP04756	FP02286	nexus 3 hub with extra drilled holes
18b	---	2996706	shifting lever
18c	---	2996707	components for hub
18d**	---	2996708	clamp 5/8"
18e	---	2996709	gear 18T (2,3mm)

* - option

** - not depicted

*** - non-modified set, non-usable for clients

2/2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Differential for 2219 with Nexus 3-speed				---
Scale:	1:7	Drafted	Verified	Confirmed
Date	2.12.2025	2.12.2025	2.12.2025	
Initials	XX	XX	XX	

Plan No.
EP00169/J
