

Bill of Material

No.	Plan	Article	Quantity
1	M8 DIN985	1996213	1
2	X002133	1992424	4
3	X003692	1998763	4
4	X004452	1904058	4
5	X004606	1905091 (39/44/48)	2
	X004607	1905092 (42/46/50)	2
6	X004620	1905090	2
7	P00023336	2639354 (39)	1
	P00023305	2639554 (42)	1
	P00023279	2638154 (44)	1
	P00023307	2639754 (46)	1
	P00023309	2639954 (48)	1
	P00023311	2640154 (50)	1
8	P00023337	2639454 (39)	1
	P00023306	2639654 (42)	1
	P00023280	2638254 (44)	1
	P00023308	2639854 (46)	1
	P00023310	2640054 (48)	1
	P00023312	2640254 (50)	1
9	3,9x9,5 DIN7981	2992349	4
10	M8 DIN126	1996214	2
11	M4 DIN985	2994901	4
12	M6 DIN985	2994063	6
13	M6x40 DIN912	2994317	4
14	M6x50 DIN912	2995258	2
15	M8x110 DIN912	2995259	1
16	M4x16	2995260	4
19	P000064	1905711	1
20	P000064	1916218	1

Sheet 1/2

V500 LIGHT CROSS

Plan No.

E00023449/F_1.1

Scale: 1:5	Drafted	Verified	Confirmed
Date	1/09/2021	13/01/2025	13/01/2025
Initials	PM	SD	PJ




Bill of Material - seat depth 42-46			
No.	Plan	Article	Quantity
17*	Seat sail W39	2997079	1
	Seat sail W42	2997080	1
	Seat sail W44	2997081	1
	Seat sail W46	2997082	1
	Seat sail W48	2997083	1
	Seat sail W50	2997084	1
18A*	B23 cushion W39	2995670	1
	B23 cushion W42	2995671	1
	B23 cushion W44	2995672	1
	B23 cushion W46	2995673	1
	B23 cushion W48	2995674	1
	B23 cushion W50	2995675	1
18B*	B24 Visco W39	2997073	1
	B24 Visco W42	2997074	1
	B24 Visco W44	2997075	1
	B24 Visco W46	2997076	1
	B24 Visco W48	2997077	1
	B24 Visco W50	2997078	1

Bill of Material - seat depth 48			
No.	Plan	Article	Quantity
17*	Seat sail W39	2995778	1
	Seat sail W42	2995779	1
	Seat sail W44	2995780	1
	Seat sail W46	2995781	1
	Seat sail W48	2995782	1
	Seat sail W50	2995783	1
18A*	B23 cushion W39	2995682	1
	B23 cushion W42	2995683	1
	B23 cushion W44	2995684	1
	B23 cushion W46	2995685	1
	B23 cushion W48	2995686	1
	B23 cushion W50	2995687	1
18B*	B24 Visco W39	1922141	1
	B24 Visco W42	1922142	1
	B24 Visco W44	1922143	1
	B24 Visco W46	1922144	1
	B24 Visco W48	1922145	1
	B24 Visco W50	1922146	1

* not depicted

Sheet 2/2

2995226

V500 LIGHT CROSS				Plan No.
				E00023449/F_1.1
Scale: 1:5	Drafted	Verified	Confirmed	
Date	1/09/2021	13/01/2025	13/01/2025	
Initials	PM	SD	PJ	