



**VERMEIREN
NV**

Vermeirenplein 1-15
2920 Kalmthout
www.vermeiren.com

Made in Poland

NeoFlex

**U01Neo.
size M**

contains:
90% Neopren
8% Poliester
2% Elastan




**SGS INSTITUT
FRESENIUS** **CE**



Material					
No.	Article	XS	S	M	L
1	2991014	0,033	0,06	0,068	0,12
2	2991886	1,5	1,7	1,8	2,2
3	1009517	1,3	1,3	1,3	1,3
4	1009512	3	3	3	3
5	1009508	2	2	2	2

Sizes [cm]			
	A	B	Article
XS	17	19	2991100
S	22	25	2991274
M	27	25	2991275
L	34	32	2991276

Mounting			
No.	Plan	Article	Quantity
1	XP01976	2993641	2
2	XP01980	2993813	2
3	XP01979	2993812	2
4	M6X30 collar bolt	2994294	2
5	YP07178	2760099	2



EP00348 / B

NEO.U01.THIGH BELT				Plan No.
				EP00348/ B
Scale: N/A	Drafted	Verified	Confirmed	
Date	31/07/2018	04/03/2019	04/03/2019	
Initials	FD	AS	TK	