



VERMEIREN

NV

Vermeirenplein 1-15
2920 Kalmthout
www.vermeiren.com

Made in Poland

NeoFlex

U80Neo.

size M

contains:

90% Neopren




8% Poliester

2% Elastan



Material						
No.	Article	XS	S	M	L	XL
1	2991014	0,011	0,013	0,015	0,018	0,023
2	1009512	2	2	2	2	2
3	1009517	0,28	0,32	0,36	0,4	0,44
4	1009791	0,35	0,4	0,45	0,5	0,55
5	1009792	0,35	0,4	0,45	0,5	0,55
6	2991886	0,7	0,8	0,9	1	1,2

Sizes [cm]		
	Ankle circ	Article
XS	14-17	2991103
S	17-20	2991102
M	19-23	2991101
L	22-29	2991277
XL	28-33	2991278



EP00350 / A

NEO U77 ANKLE STABILISATOR				Plan No.
				EP00350/ A
Scale: N/A	Drafted	Verified	Confirmed	
Date	31/07/2018	31/07/2018	31/07/2018	
Initials	FD	AS	TK	