




Bill of Material			
No.	Plan	Article	Quantity
1	Y019204	1782354	1
2	YP08073	2606254	1
3	YP08074	2606354	1
4	P00030024	2606054	2
5	Top Case	W201138	1
6	P00025186	GP27600	2
*	X005150	1788100	2
7	M8x60	1900809	4
8	-	-	-
9	M8x12	1008963	2
10	M8x25 DIN912	1990501	2
11	Washer M8	1003199	4
12	Nut M8	1001682	2
13	M6x30 DIN7984	1992721	2
14	M6x20	1006999	4
15	Washer M6	1997328	4
16	Nut M6	1001680	6
17	Tube cap Ø22	1990648	2
18	Washer M8	1001556	2

* - use with frame extender EP00497

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

TOP CASE SE27 TURIOS

Scale: 1:6	Drafted	Verified	Confirmed
Date	14/02/2020	09/06/2025	09/06/2025
Initials	PM	AZ	MK

Plan No.
EP00462/G

