



EC_COM2



No.	Plan	Article	Qty
1	P00065384	2679165	1
2	P00065219	2679265	1
3	P00065394	2679365	1
4	P00065220	2679465	1
5	-	-	-
6	XP02584	2994386	8
7	M5X20 DIN912 BLACK	1915452	4
8	WASHER M5 DIN125 A	2994069	8
9	M5x10 DIN912	2994902	4
10	M5x12 DIN912	2994124	1
11	PIN M8	2998094	1
12	M6X16 DIN7991 BLACK	2994064	2
13	-	-	-
14	M5X25 DIN912 BLACK	2994906	1
15	1996213	1996213	1

Electrical components are presented in 2998597

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

SWING AWAY LEFT PLEGO			
Scale: 1:2	Drafted	Verified	Confirmed
Date	08/05/2023	19/03/2025	19/03/2025
Initials	BK	DR	MK

2998093

Plan No.
E00066678/C.1