




No.	Plan	Article	Qty
1	Y011988	1669654	1
2	Y011987	1669854	1
3	Y011148	1669754	6
4	Lever M8x15	1006661	1
5	XP01548	1001524	2
6	Lever M8x12	1900712	3
7	Washer M8	1001554	3
8	Y006826	1599954	1
9	M6x20	1006999	2
10	X003698	1004550	2
11	Star Knob M6x20	1009932	1
12	Tube cap 15x15	1005899	1
13	Y006187	1599800	1
14	E00066476	2998020	1
15	E00052508	2998018	1
16	1001526	1001526	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

PLEGO L58C HEADREST			
Scale: 1:4	Drafted	Verified	Confirmed
Date	18.08.2023	10-04-2025	10-04-2025
Initials	BK	MP	SD

2998011

Plan No.
E00066475/C