



Material				
No.	Plan	Article	Quantity	
1	P00050415	2671226	1	
	P00050654	2671526	1	
2	P00050648	2671126	1	
	P00050656	2671426	1	
3	P00050649	2671326(39)	1	
	P00035168	2652126(42)	1	
	P00035169	2652226(44)	1	
	P00036047	2652326(46)	1	
	P00035170	2652426(48)	1	
	P00050657	2671626(50)	1	
	YP07538	2773454(39)	1	
4	YP07537	2773354(42-44)	1	
	YP07536	2773254(46-50)	1	
	* 5	P00050660	2671726(39)	1
* 5	P00035177	2653226(42)	1	
	P00035178	2653326(44)	1	
	P00035179	2653426(46)	1	
	P00035180	2653526(48)	1	
	P00050661	2671926(50)	1	
	6	M6x25 ISO7380	2994062	1
	7	M6 DIN985	2994063	1
* 8	M6 DIN125	2994068	2	
* 9	D4x14	1001584	3	

* Option

BK6 FOOTREST V500 ACTIVE				-
				Plan No.
				E00050417/A1.1
Scale:	1:3	Drafted	Verified	Confirmed
Date	2023-01-14	14/01/2023	14/01/2023	14/01/2023
Initials	AS	AZ	AS	

