




Bill of Material			
No.	Plan	Article	Quantity
1	E000810	1909057	1
2	P00087333(R)	2692354	1
*	P00087336(L)	2692454	1
3	X004415	17193	1
4	D25 Cap	1991857	1
5	M6x35 DIN7984	1903903	1
6	M6 DIN125	1001564	1
7	M6 DIN986	2994059	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

* - not depicted |
Only right side depicted

TURIOS BZ7/BZ8 FOOTREST				--
				Plan No. E00055257/B1.1
Scale: 1:3	Drafted	Verified	Confirmed	
Date	2023-03-31	17/09/2024	17/09/2024	
Initials	AS	AZ	MK	