

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

<h3>Differential for 2217 with Nexus 7-speed</h3>				---
				Plan No. EP00225/H
Scale: 1:7	Drafted	Verified	Confirmed	
Date	02.12.2016	28.11.2025	28.11.2025	
Initials	KR	PM	AZ	


Material			
No.	Plan	Article	Quantity
1	YP04707	2886600	2
2	YP04804	2992272	1
3	XP02365	GP17200	1
4	X002282	1992854	6
5	6x6x50 DIN6885A	GP19000	2
6	XP02369	GP17400	1
7	XP02368	GP17300	1
8	XP01176	2992208	1
9	XP01189	2886004	2
10	XP01202	2992224	1
11	XP03236	2631200	1
12	XP01257	2853400	1
13	XP01259	2992306 (01)	1
	XP01671	2993508 (02)	1
	XP01672	2993509 (03)	1
	XP01673	2993510 (06)	
	XP01674	2993511 (07)	1
14	XP01283	2992321	1
15	D10,5 DIN 125	1001558	3
16	M10 DIN 934	1001672	4
17	M6 DIN 985	1001680	4
18	D15 DIN 6799	1003049	4
19	EF00100 Front Wheel 20"	2991678	1
20	D160mm	2992382	1
21	D3x10 ISO 15979	2992305	6
22	M8x25	1990501	2
23	M5x16	1006112	6
24	---	1994568	6
25	M6x10	1900641	6
26	D8,4 DIN 9021	1996574	2
27	---	---	-
28	XP01274 (Nexus 7)	2992270	1
29	M10DIN 439B	1903216	1

Material			
No.	Plan	Article	Quantity
30	EF00101 Wheel 20"	2994423	2
31	XP01197	2992222	2
32**	---	1902733	3
33	Chain	2993196	1,66 [m]
34	M6 DIN 127 B	1004414	6
35	M6x90 DIN 912	1911103 (01,02,03)	4
	M6x110 DIN 912	1911104 (06)	4
	M6x100 DIN 912	1998869 (07)	4
36	XP01195	2992221	2

** - not depicted

Nexus 7 subparts - 2992270			
No.	Plan	Article	Description
28a	---	2996711	nexus 7 hub with extra drilled holes
28b	---	2996712	shifting lever
28c	---	2996713	components for hub
28d**	---	2996708	clamp 5/8"
28e	---	2996709	gear 18T (2,3mm)

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Differential for 2217 with Nexus 7-speed				---
				Plan No. EP00225/H
Scale: 1:7	Drafted	Verified	Confirmed	
Date	02.12.2016	28.11.2025	28.11.2025	
Initials	KR	PM	AZ	