




| No. | Plan           | Article | Qty |
|-----|----------------|---------|-----|
| 1   | E00066473      | 2998022 | 1   |
| 2   | Y003413        | 1463254 | 1   |
| 3   | Y007779        | 1607954 | 1   |
| 4   | Y003414        | 1463354 | 1   |
| 5   | E00052508      | 2998018 | 1   |
| 6   | Lever M8x15    | 1006661 | 2   |
| 7   | Tube cap 15x15 | 1005899 | 1   |
| 8   | M6x20          | 1006999 | 4   |

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS  
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN  
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,  
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION  
FROM VERMEIREN GROUP.

| PLEGO L55 HEADREST |            |            |            |
|--------------------|------------|------------|------------|
| Scale: 1:4         | Drafted    | Verified   | Confirmed  |
| Date               | 18.08.2023 | 17/02/2025 | 17/02/2025 |
| Initials           | BK         | WM         | MK         |

2998012

Plan No.

E00066564/B1.1

