


Material			
No.	Plan	Article	Quantity
1	YP08149	2608454 (L)	1
	YP08148	2608354 (R)	1
2	X005561	2021699	1
3	---	2997934 (L)BLACK	1
	---	2993244 (L) MOKA	1
	---	2996998 (L)Airflow	1
	---	2997935 (R)BLACK	1
	---	2993248 (R) MOKA	1
	---	2996999 (R)Airflow	1

4	XP00227	2990147	2
5	X005517	2014999	1
6	X005565	1912287	1
7	M6x45 DIN912	1003109	2
8	M5X16	1006112	2
9	---	1008459	1
10	M6x25	1006998	2
11	M5	1001679	2
12	---	2994059	3
13	---	2994069	2
14	---	3991136	2
15	M5X20 DIN912	1915452	4
16	M5 DIN985	2994006	4
17	Y011751	1662654	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Only left side depicted

SER ARMREST INOVYS 235				Plan No.
				EP00594 / B
Scale: 3:10	Drafted	Verified	Confirmed	
Date	2021-08-02	2023-06-21	2023-06-21	
Initials	AS	OY	AS	