



EP00267 / H

* - option / opcja
 ** - not depicted / nie ujeta

| | | Material | |
|------|--------------|--|----------|
| No. | Plan/Article | Description | Quantity |
| 1 | EP00247 | LOWER FRAME | 1 |
| 2* | EP00357 | CENTRAL BRAKE | 1 |
| | EP00359 | CENTRAL BRAKE FOR ONE SIDE | 1 |
| 3 | EP00256 | LEFT SIDE RAILS | 1 |
| | EP00280 | LEFT SIDE RAILS 969 | 1 |
| 4 | EP00257 | RIGHT SIDE RAILS | 1 |
| | EP00280 | RIGHT SIDE RAILS 969 | 1 |
| 5 | EP00251 | FOOT SIDE WITH ABS PLATFORM | 1 |
| | EP00276 | FOOT SIDE WITH ABS WITHOUT MOTOR | 1 |
| | EP00259 | FOOT SIDE WITH METAL LATHS | 1 |
| | EP00277 | FOOT SIDE WITH METAL LATHS WITHOUT MOTOR | 1 |
| | EP00367 | FOOT SIDE WITH METAL LATHS MOTOR 2ND POS. | 1 |
| | EP00368 | FOOT SIDE WITH ABS MOTOR 2ND POSITION | 1 |
| 6 | EP00260 | HEAD SIDE WITH ABS PLATFORM | 1 |
| | EP00261 | HEAD SIDE WITH METAL LATHS | 1 |
| 7 | EP00262 | HEAD PANEL | 1 |
| | EP00272 | HEAD PANEL 965 | 1 |
| | EP00274 | UPHOLSTERY HEAD PANEL | 1 |
| | EP00358 | HEAD PANEL FULL | 1 |
| | EP00519 | HAED SIDE PANEL / ILLICO VINTERNO | 1 |
| 8 | EP00263 | FOOT PANEL | 1 |
| | EP00272 | FOOT PANEL 965 | 1 |
| | EP00273 | UPHOLSTERY FOOT PANEL | 1 |
| | EP00358 | FOOT PANEL FULL | 1 |
| | EP00520 | FOOT SIDE PANEL / ILLICO VINTERNO | 1 |
| 9* | EP00366 | BUMPER FOR STANDARD PANEL | 1 |
| 10 | EP00268 | SELF RAISING AID | 1 |
| 11 | EP00269 | TRANSPORT SYSTEM | 1 |
| 12* | EP00386 | HORIZONTAL KIT FOR TRANSPORT SYSTEM | 1 |
| 13* | EP00371 | HORIZONTAL TRANSPORT SYSTEM | 1 |
| 14* | EP00279 | BED ELONGATION | 1 |
| 15* | EP00411 | VERTICAL TRANSPORT SYSTEM | 1 |
| 16* | EP00473 | WIDTH EXTENSION KIT | 1 |
| 17* | EP00493 | SPLIT SIDE RAILS | 1 |
| 18** | EP00278 | STICKERS | 1 |
| 19** | U000369 | CARTON FOR BED WITH STANDARD PANEL | 1 |
| | U000367 | CARTON FOR BED WITH PANEL 965 | 1 |
| | U000386 | CARTON FOR BED WITH CENTR.BR. STAND. PANEL | 1 |
| | U000387 | CARTON FOR BED WITH CENTR.BR. PANEL 965 | 1 |
| | U000424 | CARTON FOR BED WITH STANDARD PANEL VINTERNO | 1 |
| | U000425 | CARTON FOR BED WITH PANEL 965 VINTERNO | 1 |
| | U000426 | CARTON FOR BED WITH CENTR.BR. STAND. PANEL VINT. | 1 |
| | U000427 | CARTON FOR BED WITH CENTR.BR. PANEL 965 VINTERNO | 1 |
| 20** | 2993398 | CARTON CORNER | 4 |

| | | | | |
|-------------|------------|------------|------------|-----------|
| EV ILLICO | | | | Plan No. |
| | | | | EP00267/H |
| Scale: 1:20 | Drafted | Verified | Confirmed | |
| Date | 02/11/2017 | 08/06/2020 | 08/06/2020 | |
| Initials | DB | AZ | TK | |