




A	SEAT DEPTH	430, 450, 470, 490, 510
B	SEAT WIDTH	390, 450, 510
C	SEAT ANGLE	0°- 20°
D	FOOTREST LENGHT	347, 389, 424, 449, 474,
E	ARMREST HEIGHT	245, 260,275, 290, 305, 320
F	WIDTH OF ARMREST	411, 471, 531 width adjustment
G	SEAT HEIGHT	490, 520
H	FOOTPLATE HEIGHT	220
I	BACKREST ADJ	0°-30°

EP00480 /D

DIMENSIONS TURIOS ELECTRIC				Plan No.
				EP00480/ D
Scale: 1:10	Drafted	Verified	Confirmed	
Date	14/01/2020	20/03/2025	20/03/2025	
Initials	AS	AZ	MK	