




Bill of Material			
No.	Plan	Description	Quantity
1	E00054056	Side frame BI5	2
	E00057375	Side frame B02/B03	2
2	E00054636	Rear frame (44-50)	1
	E00082136	Rear frame (39-42)	1
3	E00053247	Lock mechanism	2
4	E00054057	Front fork	2
5	E00071954	Footrest B06	2
6	E00055449	Cross	1
7	E00055516	Joystick attachment	1
8	E00057691	Upholstery	1
9	E001038	Armrest BI5	2
10	2998597	Drive train standard	1
	2998566	Drive train + SE19 lights	1
11	1918834	Battery 24V 12Ah	1
	1919283	Battery 24V 17Ah	1
12*	1919284	Charger	1

* Not depicted

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

5146249

VERSO SPECIFICATION

Plan No.

E00057683/C.1

Scale: 1:8	Drafted	Verified	Confirmed
Date	11/04/2023	12/02/2025	12/02/2025
Initials	SD	DR	MK

