


PRINT IN COLOR



| #  | Reference                          | Art number | Description                                       | Qty  |
|----|------------------------------------|------------|---------------------------------------------------|------|
| 1  | EC_CON3                            | 2998579    | EC CON EW-M20-5P-3-3                              | 1    |
| 2  | EC_CON6-EC_CON8, EC_CON11-EC_CON13 | 2998612    | EC CON bootlace ferrules 1 mm x0,8 RD             | 6    |
| 3  | EC_CON5, EC_CON9, EC_CON10         | 2998610    | EC CON bootlace ferrules 0,5 mm x0,8 OR           | 3    |
| 4  | EC_FUS1                            | 2998575    | EC FUS FUSE HOLDER 5X20                           | 1    |
| 5  | EC_L1, EC_L3                       | 1914550    | EC L TAIL LIGHT BAYARD AT02L 24V                  | 2    |
| 6  | EC_L2, EC_L4                       | 1914551    | EC L FRONT LIGHT BAYARD AA02L 24V                 | 2    |
| 7  | EC_COM11                           | 1921014    | EC COM LIGHT MODULE                               | 1    |
|    |                                    | 1915936    | EC COM LIGHT MODULE SW                            | 1    |
| 8  | EC_WIRE12                          | L000028    | EC WIRE 2X0,25MM <sup>2</sup> BLACK (BROWN/WHITE) | 0,2M |
| 9  | EC_ZIP2-EC_ZIP11                   | 2993110    | EC ZIP BINDING BAND NYLON L=300MM                 | 10   |
| 10 | EC_FUS2                            | 2998574    | EC FUS FUSE 1A QUICK BLOW 5X20                    | 1    |
| 11 | BOM_EC1                            | 2998597    | BOM EC PLEGO DRIVE TRAIN                          | 1    |

BK BLACK  
RD RED  
LBL LIGHT BLUE  
BL BLUE  
LGR LIGHT GREEN  
GR GREEN  
VI VILLET  
OR ORANGE  
YE YELLOW  
WH WHITE  
LRD LIGHT RED  
PK PINK  
GY GRAY  
/ AND  
I OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS  
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN  
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,  
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION  
FROM VERMEIREN GROUP

2998582/B

**BOM EC PLEGO DRIVE TRAIN + SE19 LIGHTS**

Id: 1/1

Plan No.

2998582/B

Size: A4

Drafted

Verified

Comformed

Date:

2025-02-21

2025-02-21

2025-02-21

Initials:

FM

FM

TK

