




Bill of Material			
No.	Plan	Article	Quantity
1	P00065384	2679165	1
2	P00065219	2679265	1
3	P00065394	2679365	1
4	P00065220	2679465	1
5	-	-	-
6	XP02584	2994386	8
7	M5x20 DIN 912	1915452	4
8	M5 DIN125A	2994069	8
9	M5 DIN986	2994902	4
10	M5x10 DIN912	2994124	1
11	Pin M8	2998094	1
12	M6x16 DIN7991	2994064	2
13	Nut M8	1996213	1
14	M5x25 DIN912	2994906	1

Electrical components are presented in 2998597

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

2998023


SWING AWAY RIGHT PLEGO

Plan No.

E00065217/E.1

Scale: 1:2	Drafted	Verified	Confirmed
Date	08/05/2023	20/03/2025	20/03/2025
Initials	BK	DR	MK





NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

2998023

SWING AWAY RIGHT PLEGO

Plan No.

E00065217/E.1

Scale: 1:2

Drafted

Verified

Confirmed

Date

08/05/2023

20/03/2025

20/03/2025

Initials

BK

DR

MK

