



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Driving Gear With Nexus 3 2219			
Scale:	1:6	Drafted	Verified
Date	05.06.2014	2.12.2025	2.12.2025
Initials	AZ	PM	AZ

Plan No.
EP00156/N


Material			
No.	Plan	Article	Quantity
1	YP03014 (Fix Wheel)	2822000	1
2	-----	-----	-
3	YP04706	2886500	1
4	6x6x50 DIN6885A	GP19000	1
5	XP02362	GP17100	1
6	XP03236	2631200	1
7	XP02365	GP17200	1
8	XP01283	2992321	2
9	M10 DIN 125-1A	1001558	5
10	M10 DIN 934	1001672	6
11	D20 DIN471	1995259	4
12	5x5x40 DIN 6885A	1940800	1
13**	---	1902733	3
14	Front Wheel 24"	2992350	1
15	D160mm	2992382	1
16	EF00104 Wheel 24" Left	2994426	1
17	M5x16	1006112	6
18	M8x25	1990501	2
19	---	1994568	4
20	M8 DIN 9021	1996574	2
21	M10 DIN 439B	1903216	1
22	M5x10 DIN 912/8.8	1001488	1
23	EF00105 Wheel 24" Right	2994425	1
24	Chain	2993196	1,66 [m]
25	M4 DIN 985	1001678	6
26	XP03143	2996710	1
27	XP01177	2992207	1
28	X002282	1992854	6
29	M4x16	2990185	6
30*	---	---	-
	Vintage	1907965	3

Nexus 3 subparts - 2996710 (2992142***)			
No.	Plan	Article	Description
26a	YP04756	FP02286	nexus 3 hub with extra drilled holes
26b	---	2996706	shifting lever
26c	---	2996707	components for hub
26d**	---	2996708	clamp 5/8"
26e	---	2996709	gear 18T (2,3mm)

* - option

** - not depicted

*** - non-modified set, non-usable for clients

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Driving Gear With Nexus 3 2219			
Scale: 1:6	Drafted	Verified	Confirmed
Date	05.06.2014	2.12.2025	2.12.2025
Initials	AZ	PM	AZ

Plan No.
EP00156/N
